

# Ion Beam Profiling Using a Novel Electronic Imaging Detector

\*Bruce N Laprade and Richard A Prunier  
BURLE Electro-Optics  
Sturbridge MA

Paper 1700-2400 Presented at The 2004 Pittsburgh Conference  
March 7, 2004

# Introduction

- In the rapidly evolving analytical instrument market, new applications are constantly being developed. Twenty five years ago mass spectrometers were primarily research tools. Today, mass spectrometers are used in medical diagnostics, semiconductor manufacturing, environmental monitoring, drug discovery and food processing. Virtually everyone has benefited by the existence of these instruments.

# Introduction Cont.



- Rapidly emerging new applications for Mass Spectrometers have forced instrument manufacturers to substantially reduce product development cycle time.

# Discussion

- A number of new materials and tools including modeling software have been developed in order to facilitate rapid product development.
- Ion optics modeling software is used to design and predict the ion path within the mass spectrometer. However, ion trajectories can be influenced by many factors which are not considered in the model. Efficiently transporting ions from the source through the mass filter is critical for maximizing instrument sensitivity.

# Discussion Cont.

- The conventional method for aligning an ion beam consists of scanning the beam over a Faraday cup or electron multiplier, integrating the current, and finding the settings which produce the highest signal.
- A new imaging tool has been developed which can be used to visualize the location of any charged particle (Ion, Electron), UV photon or soft x-ray beam. Imaging the beam enables the instrument designer to ensure all available signal ions are collected.

# Typical Applications




- Ion Beam Profiling
- Ion Optic Model Verification
- Imaging TOF
- VUV Spectroscopy
- High Energy Physics

# Apparatus



- The Electronic Imaging Detector is a microchannel plate phosphor screen based detector coupled to a CCD camera with a frame grabber card and software for a personal computer.
- This device was designed to enable the user to capture images from the phosphor screen at a rate of up to 30 frames/sec and store them in a variety of formats suitable for enhancement or manipulation. TTL inputs enable event and strobe triggers to be utilized to synchronize the camera to an event.
- The internet capability of the system enables experimenters to share real time images.

# The Electronic Imaging Detector




**MCP  
Chevron™**

**6" Conflat  
Flange**

HV Feedthroughs for the Chevron and Phosphor screen



# The Electronic Imaging Detector



MCP and Phosphor  
Screen Assembly

Vacuum Flange

CCD Camera inside  
Light Shroud

# The Electronic Imaging Detector



Phosphor Screen  
High Voltage

Video Out

Chevron Out

Chevron In

Grid



Power


Vacuum Flange

CCD Camera inside  
Light Shroud

+ 9 volt in

# The Electronic Imaging Detector


CCD  
Camera




Fiber Optic  
Phosphor  
Screen

# Ion Beam Profile on a Test Pattern


---




# Chevron Resolution Test




# Random Ion Pattern



# Test Configuration for Imaging an Ion Beam From a Quadrupole Mass Filter



# Experimental Test Set-Up



Ion Source


Quadrupole Mass filter

Microchannel Plate  
Phosphor screen Detector


Vacuum Flange



# Ion Beam from a Quadrupole Mass Filter




Poorly Focused Beam



Well Focused Beam




Quad Rods in Positive Mode




Quad Rods in the Negative Mode


# Typical Reflectron Instrument Geometry




# Ion Beam Images




Beam off the detector  
Left Side



Beam off the Detector  
Top Left



Beam Centered, But not  
Focused



Loose wire found obstructing  
the beam

**PHYSICAL CHARACTERISTICS of MCPS**

Quality Diameter:	25mm Minimum
Center-to-Center Spacing:	12 $\mu$ m Nominal
Pore Size:	10 $\mu$ m Nominal
Bias Angle:	12° $\pm$ 1°
Open Area Ratio:	55% Minimum
Quality Level:	Imaging

**SPECIFICATION****ELECTRICAL CHARACTERISTICS of DETECTOR**

Electron Gain @ 2400 Volts: Maximum	1 x 10 <sup>7</sup> Minimum
Bias Current Range @ 2400 Volts:	4-29 Microamps
Resistance:	83-600 Megohms Reference
Dark Count (measured at gain voltage):	5 (cts/sec/cm <sup>2</sup> ) Maximum
Pulse Height Distribution (measured at gain voltage)	100% Maximum
Linear Output Current Density: (Microamps/cm <sup>2</sup> )	Typically 10% of Bias Current Density


**SPECIFICATION****CCD CAMERA AND FRAME GRABBER**

<b>CCD Camera:</b>	Monochrome
Pixel Size: (H X V)	7.15 X 5.55 microns
Pixels (H X V)	537 X 550
Horizontal Resolution:	380 TV Lines
Signal Format:	EIA
Features:	AGC and Auto Electronic Iris
Electronic Shutter Speed:	1/60 – 1/10,000 sec
Gamma	0.45
Minimum Sensitivity:	0.5 lux
Frame Rate:	30 FPS
Power Consumption:	0.9 watts
Image Coupling to the Phosphor Screen:	Lens
<b>Frame Grabber:</b>	
Operating System Requirements:	Windows 98, 2000, NT and XP
Hardware Installation Requirements:	Half Size PCI Slot
Installation:	True Plug and Play

**SPECIFICATIONS**

Software Installation:	CD
Bus Interface:	PCI Bus Rev. 2.1 Compliant
Max. read/write speed:	130 MBytes/sec.
Data Format:	YUV 4:2:2, RGB 565, RGB 888, XRGB and 8 bit Monochrome.
Video Input:	4 composite video, software switchable
Connector:	DB-15HD High Density Connector or RCA Connector
Frame File Type Choices:	TIFF, Bitmap and TARGA
Event Strobe and Trigger:	Two channels of TTL

The Chevron™ Model 3025FM detector assembly contains two Imaging Quality Advanced Performance Long-Life™ Microchannel Plates and a fiberoptic phosphor screen with P20 phosphor mounted to a 4.5" vacuum flange and CCD camera.



# System Features

- An ion imaging tool has been successfully developed and has been demonstrated to produce high quality images of various beam profiles.
- The microchannel plate detector provides signal amplification up to 10 million, enabling single ion or electron events to be imaged and stored.
- Automatic brightness control (ABC) and variable electronic shutter speed prevents blooming.
- The CCD camera, frame grabber and software facilitates efficient image capture, storage and manipulation.
- TTL event trigger and strobe ensures frames are collected during relevant timeframes.